



How to Install the HemaTrax[®] Product Table Patch for COVID-19 Convalescent Plasma Codes

This patch should ONLY be used for HemaTrax[®] versions v8.0 and v9.0

To install the Product Table patch:

1. Download the patch from our LSS portal using your existing account at:
<http://download.digi-trax.com>
2. Unzip the zip file.
3. Transfer all unzipped files to **C:\ProgramData\HemaTrax_Print_Server\DATA** directory to replace the existing product table files.

Please Note: If ProgramData is not listed under **C:**, unhide it. Go to: Control Panel, File Explorer, Options, click View tab, click "Show hidden files, folders, and drives" option, then click Apply and OK.

4. Go to "Control Panel - Administrative Tools" and under "Services", restart the "HemaTrax Print Server" service.

If you have any questions, please contact Henry Cagadas at 847-613-2114 or at hcagadas@digi-trax.com