



## HemaTrax CT

### 3.5.1 Release Notes

6/2019

- Product code database updated to ICCBBA version 7.25.0.
- Corrected an issue when autologous donations were appearing in the upper and lower right quadrant if the following steps occurred:
  - a) A label size other than the 4x4 was selected.
  - b) The user went to the Print Log and selected a record and pressed the Get Data Button.
  - c) On the Label Design tab, the label size was changed back to the 4x4 label
- Corrected an issue where Autologous donation GRID numbers were displaying info in lower quadrant as well.
- Corrected issue where Autologous donation Donor Date of Birth field would continuously scroll.
- Added fields Product Code Type and Expiration Type to the Print Log report.
- Adjusted reference to Daylight Savings Time to the correct term Daylight Saving Time.
- Improved alignment of the human readable text below the expiration barcode.