



HemaTrax[®]-CT 3.7.0

Release Notes

9/2021



- Product code database updated to ICCBBA version 7.56.0.
- Added support for HemaTrax[®] MFG. A new format that allows for the labeling of apheresis collection products for sponsored cellular therapy clinical trials and manufacturing, based on the ICCCBBA standard.
- Corrected a calendar scrolling bug. Calendar fields should now increment correctly when clicking the up/down arrows.
- Removed the Save Changes button in the Facilities tab. Any changes made to the facilities no longer require selecting Save Changes.
- When selecting the 48-hour expiration option with the default 23:59 time selected, the UTC time will remain at 23:59 as well. 23:59 is the ISBT 128 default time if no time is provided.
- Added “Generate DIN” activity to the activity log.
- Corrected a duplicate anonymous login that would record in the activity.