



# HemaTrax<sup>®</sup>-CT 3.7.3

## Release Notes

06/2024



### General

- Product code database updated to ICCBBA version 7.72.1.
- Clear All button now updates the label preview after clearing.
- Fixed print log search generating an error after second search with invalid value.

### HemaTrax MFG Format

- When “No Expiration” is selected, the text “No Expiration” appears on the label.
- Fixed missing “-“ separators in Donor DOB when autologous products are selected.
- “Donor DOB” and “Patient DOB” are now displayed instead of “Date of Birth”.
- Biohazard logo will now appear when selecting a biohazard collection type.
- Repositioned “Related Donor, First or Second Degree” donation type to better fit on the label.
- If Licensed Product is selected, repositioned Donor Name Donor ID and Donor DOB text so that there is no overlapping text.