





## Enhance your Process for Tracking Blood from Beginning to End

HemaTrax-BC is a new economical software offering from Digi-Trax<sup>®</sup> Corporation designed to help your facility better track the movement of blood bags in and out of the OR and ER. You know how important it is to have an easy to use, concise way of recording the transfer of blood bags back & forth from the ER, OR, refrigerated storage unit, or where ever the blood is going.




Why not simplify your current process and create a more permanent log of your facility's blood bag activity? You'll reduce reporting errors, decrease the amount of wasted blood and ensure that you are prepared for your next FDA audit. All you need is a computer and scanner that can be placed next to the storage unit. To start, scan the employee ID, select check in or check out and then scan the bag. That's all there is to it. The Blood Check Client software screen details the following:



The system automatically logs in the date scanned and the check in time. It also logs in the location of the bag and the computer. Once the blood is returned to the storage area or refrigerator, you repeat the process when the blood is next pulled. HemaTrax-BC supports both Codabar and ISBT barcode formats.

All of this information can be stored remotely on your network for future review. When you're ready to generate a report, simply go to our log viewer in HemaTrax-BC and you will see the following in real time:



For more information about the new HemaTrax-BC, contact Digi-Trax at 800-356-6126 or e-mail [info@digitrax.com](mailto:info@digitrax.com)

### About Digi-Trax -- Unparalleled Blood Labeling Expertise and Track Record

Digi-Trax has been providing blood banks and transfusion services with on-demand blood product labeling solutions since 1990. We pioneered the first FDA compliant ABO/RH and Expiration date on-demand CODABAR print solution. Our staff consulted on the ISBT-128 Task Force that developed the U.S. Consensus standard. We remain a continuing member of America's Technical Advisory Group (ATAG) that manages the ISBT specification. Our labeling systems are in use at most of the nation's major blood banks as well as all U.S. military blood donation centers. Digi-Trax is a registered ICCBBA software developer and AABB member.

