



Cord Blood Product Labels

It's a once in a lifetime process for parents to store their child's blood

The solution for your cord blood labeling

We understand how important it is to properly identify and manage the cord blood that you receive.

Digi-Trax® has been supplying blood banks, hospitals and healthcare facilities with durable, smudgeproof labels and labeling systems to properly identify cord blood unit bags and testing sections. Our labels withstand extreme temperatures, both high and low, even during a vapor phase and liquid nitrogen process. And we understand that these labels may need to be stored indefinitely so our R&D team does extensive testing to ensure the labels work in these environments.



Unprocessed cord blood

Cord Blood Segment Label
product # **CRYO-CBA-CG-KT**
Size: 2.5" x 3.125"

For more details, go to:
www.digi-trax.com/ct-labels

Take a look at what one cord bank facility does

Label sample after it's printed and in use:



Labels for segments 1, 2 and 3 are folded along the tube

For more information about cord blood labeling, contact Digi-Trax at 800-356-6126 or email us at info@digi-trax.com.